CREA

Hersteller

Schweiz

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Machine read, semi-quantitative detection of creatinine in urine

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Package contents:

10 slide tests, singularly welded in aluminium foil each with desiccant, overflow chamber pipette and work instructions.

Reagents:

Supplied in addition with the test:

- Creatinine dilution buffer:
- 900 µl buffer for deluting 100 µl urine sample
- Barcode for P.I.A.²

Additionally required materials:

- Urine cup
- Timer
- Gloves

Optionally available materials:

- Urine cup with screw lid
- Gloves
- Disposal bag
- Thermal strips for temperature control
- Result protocols for documentation

Intended use:

The Crea quick test is a sensitive immunoassay for semi-quantitative detection of creatinine in urine with P.I.A.² Systems of Protzek. It serves for detecting internal or external manipulation of urine by diluting or replacing urine with another solution.

Introduction:

Normally, a urine check includes a pre-screening of urine via immunoassay to proof drug abuse and when the result is positive, an evidencing confirmation analysis of the corresponding sample with mass spectrometry (in road traffic also a blood sample, if applicable). A negative pre-test that would have to be confirmed accordingly, is often omitted, though. To avoid confirmation analysis, test persons often try to manipulate the sample such that a false negative result is achieved and no confirmation analysis is performed, for which evidence could be provided despite manipulation.

Possible manipulations are:

- Externally diluting already discharged urine with water or with yellow tinted aqueous fluids.
- Internally diluting the drug amount in urine through increased urination due to excessive drinking or consumption of diuretics (substances to stimulate urine production).
- Disturbing the analysis principle with substances added to the urine (e.g. vinegar-based cleaner, baking soda, drain cleaner etc.).

Creatinine is a decomposition product of creatine in the muscle. It is normally created in the body (depending on the muscle mass) at a quite constant rate and is excreted via the kidneys. Manipulations by diluting urine, increased urination or replacing urine with other fluids can be detected with the Crea test through lowered Creatinine levels. Under normal circumstances the creatinine levels are between 80 and 200 mg/dll.^{1, 2, 3, 4, 5} Urine samples with creatinine levels of under 20 mg/dl are considered as not usable.

Summary:

The Crea test serves as a tool to proof manipulated urine by diluting or excessive drinking or by replacing urine with another fluid.

Test principle:

The Crea test is based on a specific antigen antibody reaction and is a competitive immunochromatographic test.

Each test contains the colour-marked specific antibody against creatinine (gold linked, monoclonal mouse antibody), the membrane-bound control antigen (Crea-BSA) and a second membrane-bound antibody for test control (polyclonal anti-mouse IgG antibody). The creatinine in the sample competes with the creatinine molecules fixed on the test membrane for a limited number of specific binding sites at colour-marked antibodies. Accordingly, the test band will be brighter the more creatinine is contained within the sample. There is no test line, when the antibodies are saturated with creatinine.

The Crea test can be evaluated with the analytical instruments of the P.I.A. product line. Formation of an antigen antibody gold complex at the test line is detected and quantified by the device. The device additionally provides an evaluation of the numeric result.

Important notes

1. Use package content only for in vitro diagnostics.
2. Keep away from children.
3. Do not swallow and avoid direct contact of sample materials and test reagents with skin and mucous membranes.
4. Urine samples are considered potentially infectious materials. Please observe standard hygienic guidelines when dealing with body fluids and have disinfectants ready.
5. The performed test can be disposed using the disposal bag. Please comply with national and regional regulations.
6. Do not use test after the expiry date printed on it.

Sample collection:

Collect urine samples in plastic or glass containers. Use fresh samples (sample is warm after discharge). Samples should be collected under visual inspection to avoid manipulation.

Storage of tests and samples:

Store test only in the original packaging and at **temperatures between 35°F and 46°F; +2°C and +8°C (in a refrigerator)**. Avoid direct sunlight. Do not freeze tests! Samples can be preserved up to 48 hours in the refrigerator (35°F to 46°F; +2°C bis +8°C). Freezing samples (-4°F; -20°C) is possible. After defrosting bring samples to room temperature and mix well.

Test execution:

All reagents are included in the test ready for use. The enclosed buffers are required to execute the test. Test handling temperature: 59°F to 86°F (+15°C to 30°C).

1. Take required test out of the packaging immediately prior to executing the test.
2. Document the test person's identity (ID) on the test.
3. Optional: Stick thermal strip in the area of the urine fill level on the urine cup and read. Colour shows the temperature in °C. This should be in the range of 32.5°C to 37.7°C (90.5°F to 99.9°F).⁶
4. Enter/select test person ID and user in the device and insert barcode for your Crea test with the barcode facing down into drawer.
5. Open a buffer tube and have it ready.
6. Take the pipette out of the packaging and pull up urine from the cup. It is important that the urine flows into the overflow chambers.
7. Add urine into the buffer tube, close it tightly and mix well.
8. Pull up the diluted sample with the pipette. It is important that the sample flows into the overflow chambers. Apply the entire pipette content onto the test.
9. Start timer and start measurement exactly after 15 minutes with the check button (V) at the device.

Alternatively: Start the set timer with the button "8".

What the results reveal:

mg/dl	Category	Interpretation
< 20	Urine sample not usable	Creatinine concentration is extremely low. Request another sample possibly after 1 to 2 hours and evaluate.
20 – 100	Low	Low creatinine content, urine diluted through high fluid addition → measured drug concentration accordingly higher.
100 – 200	Normal	Measured drug concentration corresponds approximately with actual amount.
200 - 400	High	High creatinine content, urine concentrated through too little fluid addition → measured drug concentration accordingly lower.
400 - 1000	Very high	Creatinine concentration is very high → possibly falsified through taking creatinine → measured drug concentration accordingly lower.

Restrictions and interferences:

The Crea quick test is only to detect improper dilution i.e. designed to provide information about the discharge status. Other changes or manipulations are not detected. The Crea test is not applicable to diagnose a kidney insufficiency and other illnesses. There are currently no reports available regarding possible cross reactions. Regarding the latter we refer to our constantly updated medical scientific information that can be obtained from us.

Limits and ranges:

Measuring range:

10 to 1000 mg/dl

Sensitivity:

10 mg/dl.

Sensitivity (lower detection limit) is defined as lowest analyte concentration that still can be differentiated by the analyte free solution. It is calculated as the concentration that is three standard deviations under the lowest standard (n=10).

Literature

1. Edwards C, Fyfe MJ, Liu RH and Walia AS. *Evaluation of Common Urine Specimen Adulteration Indicators*. J Anal Toxicol 17, 241 (1993)
2. Murray RL. *Creatinine*. In *Clinical Chemistry: Theory, Analysis and Correlation*, 2nd Edition, Kaplan LA and Pesce AJ (Eds.) Pages 1015-1020. (1989)
3. Newkirk RE and Rawnsley HM, *Creatinine Clearance*, SCP Check Sample Clinical Chemistry, No. CC-110, Chicago, American Society of Clinical Pathologists (1978)
4. Faulkner WR and King JW, *Renal Function*. In Tietz NW, Ed. *Fundamentals of Clinical Chemistry*, No. CC-110, Chicago, American Society of Clinical Pathologists (1978)
5. Needleman S.B., Porvaznik M., and Ander D. *Creatinine Analysis in Single Collection Urine Specimens*. Journal of Forensic Sciences 37, Seiten 1125-1133 (1992)
6. *Mandatory Guidelines for Federal Workplace Drug Testing Programs*. National Institute on Drug Abuse. Federal Register Vol. 53, No. 29, 11979 (1988)